LAICA®

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Type VP-M7

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Prodotto da: Vapo Healthcare Co., Ltd
Southern Unit of Third Floor, Building B, No.99 Yudai West Rd,
High Tech District, Kunshan, Suzhou, 215301, Jiangsu, P.R. China
Share Info GmbH
Am Schulzentrum 12, 41564 Kaarst, Germany

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STANDARDS

Il prodotto risponde ai seguenti Standards: 93/42/EEC - 1993+Amd:2007 - Medical devices directives, EU 2017/745, EN ISO 14971 - 2019 - Medical devices-Application of risk management to Medical devices, EN ISO 13485 - 2016+A1:2021 - Medical devices - Quality management systems - Requirements for regulatory purposes, ISO 15223-1 - 2021 - Medical devices – Symbols to be used with information to be supplied by the manufacturer Part1:General Requirement, EN ISO 20417 - 2021 - Medical devices. Information to be supplied by the manufacturer, EN 60601-1 - 2006/A2:2021 - Medical Electrical Equipment - Part 1: General Requirements for Basic Safety and Essential Performance, EN 60601-1-2 - 2015/A1:2021 - Medical Electrical Equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard:Electromagnetic Compatibility - Requirements and tests, EN 60601-1-6 - 2010+A1:2015+A2:2021 - Medical electrical equipment- Part 1-6: General requirements for basic safety and essential Performance - Collateral standard: Usability, EN 60601-1-11 - 2015+AD1:2020 - Medical Electrical Equipment - Part 1-11: General requirements for basic safety and essential performance - Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment, EN 62366 - 2015/A1:2020 - Medical devices- Application of usability engineering to medical devices, MEDDEV 2.7.1 - 2016 - GUIDELINES ON MEDICAL DEVICES - CLINICAL EVALUATION-A GUIDE FOR MANUFACTURERS AND NOTIFIED, ISO 27427 - 2023 - Anaesthetic and respiratory equipment - Nebulizing systems and components, EN ISO 10993-1 - 2020 - Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process, EN ISO 10993-5 - 2009 - Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity, EN ISO 10993-10 - 2023 - Biological evaluation of medical devices - Part 10: Test for irritation and skin sensitization, EN ISO 10993-23 - 2021 - Biological evaluation of medical devices - Part 23: Tests for irritation, ISO 18562-2 - 2017 - Biocompatibility evaluation of breathing gas pathways in healthcare applications – Part 2: Tests for emissions of particulate matter, ISO 18562-3 - 2017 - Biocompatibility evaluation of breathing gas pathways in healthcare applications – Part 3: Tests for emissions of volatile organic compounds (VOCs), ISO 18562-4 - 2017 - Biocompatibility evaluation of breathing gas pathways in healthcare applications – Part 4: Tests for leachables in condensates.

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Working principles:

Device for ultrasonic aerosol therapy with a very refined design and new technology, characterized by almost silent operation. With the traditional ultrasonic aerosol devices, compounds for inhalation are transformed, by means of oscillation of a piezoelectric crystal embedded in the device, into liquid micro-drops dispersed in air. In this aerosol therapy device, this transformation occurs by means of the flow of the drug in the micro-holes of a foil. The microdroplets of the drug, approximate size 5 µm, are able to penetrate deep into the upper and lower respiratory tract, targeting the active ingredients where needed.

Aerosol therapy is effective for the treatment of lung disorders, tracheitis, asthma, bronchitis, sinusitis and allergies. It can also instantly relieve symptoms caused by colds and allergic irritations, decreasing the urge to cough, and alleviating sinus pain and pressure, facilitating expectoration. It is an excellent device for home use, is practical and easy to use. It is built in compliance with current European standards regarding general requirements for the safety of medical electrical equipment (Standard EN 60601-1).

The product's small size makes it extremely convenient for those who must use it several times a day and for those who travel. The use of oil-based drugs or those with solid elements as active ingredients (suspensions), viscous solutions or essential oils, is not recommended. **Always follow the doctor's instructions for the type of medicine to use, the dosage, frequency and duration of inhalations.**

BEFORE USE IT IS IMPORTANT TO READ THE INSTRUCTIONS AND WARNINGS CONTAINED IN THIS BOOKLET CAREFULLY AND KEEP IT IN A SAFE PLACE FOR FUTURE REFERENCE.

KEY TO SYMBOLS

CE0197

In compliance with European legislation on medical devices

Manufacturer

Symbol for "Type BF applied parts"

LOT

Batch production number

EC

REP

European Representation

Class II

Caution! Read the user instructions carefully

Warning!

SN

Serial number

MD

Indicates the given product is a medical device

Production date

IPX2: Degree of protection of electrical equipment covers, where the first figure indicates the degree of protection against penetration by solid foreign bodies (from 0 to 6) and the second figure the degree of protection against penetration by liquids (from 0 to 8).

SAFETY WARNINGS

Using any electrical appliance entails compliance with a number of fundamental rules:

- Unplug the appliance from the mains straight after use and before commencing cleaning or inserting/removing accessories.
- Do not unplug by tugging on either power cord or appliance.
- If an electrical appliance falls into the water, do not attempt to retrieve it; instead, unplug immediately.
- Do not use the appliance while bathing or showering.
- Never immerse the appliance or USB cable in water, hot water, or other liquids. This may cause electrical leakage, electric shock, or device malfunction.
- Never leave the machine running unattended. Do not allow the appliance to be used by children, invalids or disabled individuals without adult supervision.
- This appliance must be used solely for the purpose for which it was designed and in accordance with the operating instructions. Any other use shall be considered improper and hence hazardous. The manufacturer cannot be held responsible for any damage resulting from improper or incorrect use.
- Keep all parts of the packaging away from children (choking hazard).
- Kids and infants may swallow small pieces. The main unit may cause injury to the children or infants. In case a child swallows small pieces, immediately consult your doctor.
- Before using the appliance, make sure it is intact and does not feature any visible signs of damage. If in any doubt, do not use the appliance and contact your dealer.
- In the event of a fault and/or irregular operation, turn the appliance off at once without tampering with it. Always contact your dealer for repairs.
- Keep the appliance away from sources of heat.
- Make sure your hands are dry when unplugging or plugging in the appliance and when turning the function control knob.
- This device shall be used only with the voltage stated in this booklet.
- Device with enclosure not protected against the penetration of liquids.
- Do not drop medication on the main unit. If you drop medication, immediately wipe it off with gauze. Using the unit whilst wet it may cause damage or result in electric shock.
- Do not poke the metal mesh with a cotton swab or pin, as this may damage the unit and render it unusable.
- Use only the supplied USB cable, do not use a damaged cable.
- Do not disassemble, repair, or modify the unit.
- Do not drop or expose the main unit or nebulization chamber to strong impacts or shocks. This may cause damage, rendering the unit unusable, or result in electric shock.
- Do not use the device close to high-frequency electromagnetic emitters.
- Do not use the appliance in the presence of flammable gases and in very dusty areas.

- Do not use the device, if in the atomisation chamber there is not liquid or if the liquid is present in limited quantities.
- Do not use powder medicinal products, powder medicinal products associated to liquid ones, or acid medicinal products with this appliance. Do not atomize water, as it may worsen your symptoms. Always use saline.
- Always clean the appliance after each use in order not to let the residual medicinal products solidify and damage the mechanism.
- During inhalation, keep the unit away from the eyes because some medicines can irritate and / or burn eyes.
- The device must only be used with the original accessories indicated in the instruction manual.
- Do not use the device in a room where a spray was previously used. If you still wish to use the unit, we recommend airing the room before starting treatment.

⚠ IMPORTANT WARNINGS

- In case of diabetes or other illnesses, it is suggested to consult your physician before undergoing the inhalation therapy.
- Pregnant women should consult their physician before undergoing the inhalation therapy.
- Before using a preparation, clear with your physician, if it is suitable for the person to use it. Use only medicines prescribed or suggested by your physician.
- To dilute the preparation, use exclusively an isotonic solution of sodium chloride, since otherwise respiratory problem may arise.
- This unit is designed for personal use only and should not be shared with others.
- The nebulization characteristics of this unit vary depending on the properties of the medication. Medications with high surface activity or viscosity, such as solubilizing agents or expectorants, may reduce the atomisation rate. Additionally, low medication temperatures may also decrease the nebulization rate.

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TECHNICAL SPECIFICATIONS

- Product name: Ultrasonic nebuliser for aerosol therapy
- Commercial name: NE1012
- Power source: built-in rechargeable 3.7V 800mAh lithium battery
- Power consumption: approximately 2W
- Power supply: USB cable 5V/1A, type C
- Battery life: can be recharged over 300 times
- Nebulisation chamber life: 180 hours
- Nebulization rate: ≥ 0.20 ml / min
- Particle size: MMAD approximately 5 μ m
- Capacity of atomisation chamber: ≤ 8 ml
- Pollution: Degree 2
- Voltage Category: CATEGORY II

Operating conditions:

- Temperature: min 10°C; max 40°C
- Humidity: 30%-85% RH
- Altitude: ≤ 2000 m

Storage conditions:

- Temperature: min. -20°C; max. 50°C
- Humidity: 30%-85% RH
- Atmospheric Pressure: 86kPa~106kPa

INSTRUCTIONS FOR USE

- 1) **Before each use, thoroughly clean and disinfect** the atomisation chamber, lid, dispenser, electrodes, mouthpiece and masks according to the paragraph "Cleaning and Maintenance".
- 2) Place the device on a flat surface.
- 3) Place the atomisation chamber on the main unit and it slide back until a "click" is heard (fig 1).
- 4) Fill the atomisation chamber:
 - Open the atomisation chamber lid by pulling the small hook slightly forwards and the lid upwards.
 - Add the medication, ensuring not to exceed the maximum 8ml capacity (fig 2).

- Close the lid firmly, ensuring it clicks into place.

Follow the product's instructions for use and dosage. If the product needs to be diluted with an isotonic sodium chloride solution, first add the solution, then the medication, ensuring the total volume does not exceed the maximum capacity. Exceeding the capacity may cause liquid overflow and damage the device.

Caution! The device can only be used if medication is poured into the atomisation chamber!

5) Select the appropriate attachment:

- For lower respiratory tract disorders (e.g. cough or bronchial infections), use the mouthpiece.
- For nasal and pharyngeal treatments, or for children under 5 years old, use the mask.

6) Attach the mouthpiece or adult/child mask to the dispenser (fig. 3).

7) Sit upright and hold the unit in your hand.

If using the mouthpiece, hold it firmly between the lips.

If using the mask, gently press it against your nose and ensure it covers both the mouth and nose without being too tight.

8) Turn on the device by pressing the "O/I" button. The indicator light will flash green during nebulization. Inhale slowly and deeply through your mouth, then exhale through your nose. To enhance therapy, hold your breath for a few seconds after inhaling. Remain calm and relaxed during the treatment.

If you need a brief break, remove the mouthpiece or mask and reapply it when ready to continue.

Caution: Do not shake the device when the atomisation chamber is full to avoid damage.

9) Turn off the unit using the "O/I" button when the atomisation flow has stopped.

10) Allow the unit to cool for a few minutes and proceed with cleaning as indicated in the paragraph "Cleaning and maintenance" after each treatment.

Caution! Medication may accumulate on the internal surface of the mouthpiece during usage.

If this occurs, follow these steps:

- Disconnect the mouthpiece from the dispenser carefully to avoid spilling any remaining medicine.
- Pour any medication back into the atomisation chamber, or wipe small amounts off with clean gauze.
- Inhale any remaining medication (if more than 0.5ml).

Important Notes:

- Inhalation sessions should not exceed 15 minutes. If the treatment isn't completed, press the "O/I" button to stop.
- To prolong the life of the atomisation mechanism, allow the unit to cool for an hour before continuing the treatment if it is paused.
- It's normal for some medication to remain in the atomisation chamber after the treatment, referred to as the "residual volume." This cannot be atomised.
- The device has an automatic shut-off feature when the medication is nearly used up. Never attempt to operate the unit when empty, as this can cause damage.

Attention:

- The unit requires recharging when the indicator light flashes orange. Use the USB cable supplied to recharge the unit. The indicator light will change from orange to green when fully charged.
- If the charge is low, the unit will turn off automatically after a short period of time, even if there is still some inhalation solution in the nebulisation chamber.
- Recharging will stop if the "O/I" button is used to switch on the device.

The atomisation chamber is made from consumable material, and its useful life depends on factors like usage time, exposure to oily medications, and how well cleaning and maintenance are followed.

The atomisation chamber is designed to last for approximately one year with 30 minutes of daily use.

Microbial Contamination Warning:

This unit is intended for personal use only. It is not suitable for use by multiple users.

CLEANING

1) Attention! Hot surfaces! Before cleaning the unit, allow it to cool as the transducer and residual inhalation preparations become hot during use.

Always clean the unit after each use so that the residual medicine does not solidify and ruin the mechanism.

- 2) Before cleaning switch off the device and, if recharging, unplug.
- 3) After each use, remove the atomisation chamber from the main unit by pressing the release button (fig 4). Detach the mouthpiece or mask from the dispenser.
- 4) Lift the lid on the atomisation chamber to empty any residues from the chamber.
- 5) Wash inhalation mask or mouthpiece, atomisation chamber in water. Carefully rinse the metal mesh when using medication with a high surface activity or viscosity, such as medication solubilizing agent or expectorants, as they tend to leave residues on the mesh. Do not to touch the metal mesh directly to avoid damage. Mesh cannot be used once damage occurs.
- 6) Immediately after washing the parts, wipe off water with new gauze and dry them thoroughly.
Caution! Never try to wipe off water with tissues or cloth. Dust or lint left on the metal mesh may affect nebulization.
- 7) Do not immerse the main unit in water. Instead, wipe the outside of the unit with a soft, damp cloth and mild soap. Make sure to keep the electrodes clean and dry.
- 8) Before hooking the atomisation chamber on the main unit, make sure it is perfectly dry.
- 9) We recommend disinfecting the accessories after each use to limit the risk of infections.
 - Disinfection of atomisation chamber: Soak the atomisation chamber in 75% medical alcohol for at least 3 minutes, then dry in a clean environment. The atomization chamber should be carefully stored after it is dried. Never try to wipe off water with tissues or cloth. Dust or lint left on the metal mesh may affect nebulization. Do not poke the metal mesh with a cotton swab or a pin. Do not to touch the metal

mesh directly to avoid damage. Mesh cannot be used once damage occurs. Dry it naturally.

- Disinfection of mouthpiece and mask: thoroughly clean the surface of the mouthpiece and mask with 75% medical alcohol. They should be carefully stored in a clean environment after they are dried.

Disinfect all of the components before using the appliance for the first time.

- 10) The appliance, atomisation chamber and all accessories should be perfectly dry before use.
- 11) Store the appliance and all accessories in the provided case, in a cool dry location.

TROUBLESHOOTING

Problem	Possible cause	How to correct
The nebulization rate is extremely low.	The batteries are low/worn out.	Recharge the batteries.
	The Metal Mesh is stained.	Clean the chamber to remove the stain. If this does not resolve the issue, replace the chamber with a new one.
The unit does not nebulize.	The batteries are low.	Recharge the batteries.
	The USB Cable is not correctly connected.	Remove the USB Cable from the main unit and reconnect it correctly.
	The chamber is not installed correctly to the main unit.	Install the chamber correctly.
The backlight of O/I button works correctly but the unit does not nebulize.	The Metal Mesh is stained severely.	Replace the chamber with a new one.
	The Metal Mesh is broken.	Replace the chamber with a new one.
	Medication or water is collected on the electrodes.	Clean the electrodes.
	The electrode is stained.	Clean the electrodes.

If the unit does not nebulize normally after taking the above-mentioned procedure, contact the store where you purchased the unit or the nearest dealer.

DISPOSAL PROCEDURE (Dir. 2012/19/Eu-WEEE)

The symbol on the bottom of the device indicates the separated collection of electric and electronic equipment. At the end of life of the device, do not dispose it as mixed solid municipal waste, but dispose it referring to a specific collection centre located in your area or returning it to the

distributor, when buying a new device of the same type to be used with the same functions. If the appliance to be disposed of is less than 25 cm, it can be returned to a retail location that is over 400 m² without having to purchase a new, similar device.

This procedure of separated collection of electric and electronic devices is carried out forecasting a European environmental policy aiming at safeguarding, protecting and improving environment quality, as well as avoiding potential effects on human health due to the presence of hazardous substances in such equipment or to an improper use of the same or of parts of the same.

Caution! The wrong disposal of electric and electronic equipment may involve sanctions.

WARRANTY

This appliance is guaranteed for 2 years from the time of delivery of the goods, or another longer term envisaged by the national legislation of the consumer's residence. This provision complies with Italian and European legislation. The Laica products are designed for home use and must not be used in public venues. The warranty only covers manufacturing defects and does not apply if the damage is caused by an accidental event, incorrect use, negligence or misuse of the product. Use only the accessories supplied; the use of different accessories may result in invalidity of the warranty. Do not open the unit for any reason; in the case of opening or tampering, the warranty is definitively voided. This warranty does not apply to parts subject to wear or to the batteries when supplied. After 2 years from delivery, or another longer term envisaged by the national legislation of the consumer's residence, the warranty expires; in this case, the technical assistance interventions will be carried out against a fee. Information on technical assistance, whether under warranty or for a fee, can be requested by contacting info@laica.com. No payment will be due for repairs or replacements of products that fall within the terms of the warranty. In the event of faults, contact the retailer. Do NOT send the appliance directly to LAICA. All the operations under warranty (including those of replacement of the product or part thereof) will not prolong the duration of the original period of warranty of the product replaced. The manufacturer declines any liability for any damage that may, directly or indirectly, be caused to persons, property or animals as a result of the non-observance of all the requirements established in the relevant instructions manual and concerning, especially, warnings relating to installation, use and maintenance of the appliance. Laica, in its constant commitment to improving its products, is entitled to changing without any notice, in whole or in part, its products in relation to production requirements, without this entailing any liability for Laica towards its dealers. For further information: www.laica.it.

Produced by Vapo Healthcare Co., Ltd

Southern Unit of Third Floor, Building B, No.99 Yudai West Rd,
High Tech District, Kunshan, Suzhou, 215301, Jiangsu, P.R. China

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Phone: +39.0444.795314 info@laica.com - www.laica.it - Made in China

STANDARDS

The product corresponds to the below standards: 93/42/EEC - 1993+Amd:2007 - Medical devices directives, EU 2017/745, EN ISO 14971 - 2019 - Medical devices-Application of risk management to Medical devices, EN ISO 13485 - 2016+A1:2021 - Medical devices - Quality management systems - Requirements for regulatory purposes, ISO 15223-1 - 2021 - Medical devices — Symbols to be used with information to be supplied by the manufacturer Part1:General Requirement, EN ISO 20417 - 2021 - Medical devices. Information to be supplied by the manufacturer, EN 60601-1 - 2006/A2:2021 - Medical Electrical Equipment - Part 1: General Requirements for Basic Safety and Essential Performance, EN 60601-1-2 - 2015/A1:2021 - Medical Electrical Equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic Compatibility - Requirements and tests, EN 60601-1-6 - 2010+A1:2015+A2:2021 - Medical electrical equipment- Part 1-6: General requirements for basic safety and essential Performance - Collateral standard: Usability, EN 60601-1-11 - 2015+AD1:2020 - Medical Electrical Equipment - Part 1-11: General requirements for basic safety and essential performance - Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment, EN 62366 - 2015/A1:2020 - Medical devices - Application of usability engineering to medical devices, MEDDEV 2.7.1 - 2016 - GUIDELINES ON MEDICAL DEVICES - CLINICAL EVALUATION-A GUIDE FOR MANUFACTURERS AND NOTIFIED, ISO 27427 - 2023 - Anaesthetic and respiratory equipment - Nebulizing systems and components, EN ISO 10993-1 - 2020 - Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process, EN ISO 10993-5 - 2009 - Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity, EN ISO 10993-10 - 2023 - Biological evaluation of medical devices - Part 10: Test for irritation and skin sensitization, EN ISO 10993-23 - 2021 - Biological evaluation of medical devices - Part 23: Tests for irritation, ISO 18562-2 - 2017 - Biocompatibility evaluation of breathing gas pathways in healthcare applications — Part 2: Tests for emissions of particulate matter, ISO 18562-3 - 2017 - Biocompatibility evaluation of breathing gas pathways in healthcare applications — Part 3: Tests for emissions of volatile organic compounds (VOCs), ISO 18562-4 - 2017 - Biocompatibility evaluation of breathing gas pathways in healthcare applications — Part 4: Tests for leachables in condensate.

Guidance and Manufacturer's Declaration

- 1) This equipment needs to be installed and put into service in accordance with the information provided in the ACCOMPANYING DOCUMENTS;
- 2) IMMUNITY TEST LEVELS for basic safety and essential performance of ME equipment and ME systems should be chosen based on a high probability of maintaining basic safety and essential performance, and shall be according to the professional Healthcare facility environment, home Healthcare environment, and special environment, based on the locations of intended use.
- 3) HOME HEALTHCARE ENVIRONMENT is dwelling place in which a patient lives or other places where patients are present, excluding professional Healthcare facility environments where operators with medical training are continually available when patients are present. Such as schools, outdoors, domiciles, vehicles hotels and pensions.

EXAMPLE: As indicated in Table 6 of IEC 60601-1-2: 2014 for ME EQUIPMENT, at typical cell phone with a maximum output power of 2W yields $d = 3.3m$ at an IMMUNITY LEVEL of 3V/m.

A1 Electromagnetic Emissions-For all EQUIPMENT and SYSTEMS

Guidance and manufacture's declaration-electromagnetic emission		
The DEVICE is intended for use in the electromagnetic environment specified below. The customer or the user of the DEVICE should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The DEVICE uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 14	Class B Conforms	It is possible to use the device in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings user for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A Conforms	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

A2 Electromagnetic Immunity -For Home Healthcare Environment EQUIPMENT and SYSTEMS

Guidance and manufacture's declaration – electromagnetic immunity			
The DEVICE is intended for use in the electromagnetic environment specified below. The customer or the user of the DEVICE should assure that it is used in such an environment.			
Emissions test	Test level EN 60601-1-2	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) EN 61000-4-2	±8kV contact ±2kV, ±4kV, ±8kV, ±15kV air	±8kV contact ±2kV, ±4kV, ±8kV, ±15kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Burst/Fast Transient EN 61000-4-4	±2kV power supply lines	±2kV power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge EN 61000-4-5	±1kV differential mode ±2kV common mode	±1kV differential mode ±2kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines EN 61000-4-11	0 % UT; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°	0 % UT; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions. It is recommended that the DEVICE be powered from an uninterruptible power supply or a battery.
Power frequency magnetic field EN 61000-4-8	30A/m	30A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Conducted RF IEC 61000-4-6	3Vrms 150kHz to 80MHz	3Vrms	The MANUFACTURER should consider reducing the minimum separation distance, based on RISK MANAGEMENT, and using higher IMMUNITY TEST LEVELS that are appropriate for the reduced minimum separation distance. Minimum separation distances for higher IMMUNITY TEST LEVELS shall be calculated using the following equation: $E = \frac{6}{d} \sqrt{P}$
Radiated RF IEC 61000-4-3	10V/m 80MHz to 2.7GHz	10V/m	Where P is the maximum power in W, d is the minimum separation distance in m, and E is the IMMUNITY TEST LEVEL in V/m.

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Producido por: Vapo Healthcare Co., Ltd
Southern Unit of Third Floor, Building B, No.99 Yudai West Rd,
High Tech District, Kunshan, Suzhou, 215301, Jiangsu, P.R. China

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NORMAS

El producto cumple las siguientes normas 93/42/EEC - 1993+Amd:2007 - Medical devices directives, EU 2017/745, EN ISO 14971 - 2019 - Medical devices-Application of risk management to Medical devices, EN ISO 13485 - 2016+A1:2021 - Medical devices - Quality management systems - Requirements for regulatory purposes, ISO 15223-1 - 2021 - Medical devices - Symbols to be used with information to be supplied by the manufacture Part1:General Requirement, EN ISO 20417 - 2021 - Medical devices. Information to be supplied by the manufacturer, EN 60601-1 - 2006/A2:2021 - Medical Electrical Equipment - Part 1: General Requirements for Basic Safety and Essential Performance, EN 60601-1-2 - 2015/A1:2021 - Medical Electrical Equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard:Electromagnetic Compatibility - Requirements and tests, EN 60601-1-6 - 2010+A1:2015+A2:2021 - Medical electrical equipment - Part 1-6: General requirements for basic safety and essential Performance - Collateral standard: Usability, EN 60601-1-11 - 2015+AD1:2020 - Medical Electrical Equipment - Part 1-11:General requirements for basic safety and essential performance - Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment, EN 62366 - 2015/A1:2020 - Medical devices - Application of usability engineering to medical devices, MEDDEV 2.7.1 - 2016 - GUIDELINES ON MEDICAL DEVICES - CLINICAL EVALUATION-A GUIDE FOR MANUFACTURERS AND NOTIFIED, ISO 27427 - 2023 - Anaesthetic and respiratory equipment - Nebulizing systems and components, EN ISO 10993-1 - 2020 - Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process, EN ISO 10993-5 - 2009 - Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity, EN ISO 10993-10 - 2023 - Biological evaluation of medical devices - Part 10: Test for irritation and skin sensitization, EN ISO 10993-23 - 2021 - Biological evaluation of medical devices - Part 23: Tests for irritation, ISO 18562-2 - 2017 - Biocompatibility evaluation of breathing gas pathways in healthcare applications - Part 2: Tests for emissions of particulate matter, ISO 18562-3 - 2017 - Biocompatibility evaluation of breathing gas pathways in healthcare applications - Part 3: Tests for emissions of volatile organic compounds (VOCs), ISO 18562-4 - 2017 - Biocompatibility evaluation of breathing gas pathways in healthcare applications - Part 4: Tests for leachables in condensate.

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STANDARDS

Das Produkt entspricht den folgenden Normen: 93/42/EEC - 1993+Amd/2007 - Medical devices directives, EU 2017/745, EN ISO 14971 - 2019 - Medical devices-Application of risk management to Medical devices, EN ISO 13485 - 2016+A1:2021 - Medical devices - Quality management systems - Requirements for regulatory purposes, ISO 15223-1 - 2021 - Medical devices - Symbols to be used with information to be supplied by the manufacturer Part1:General Requirement, EN ISO 20417 - 2021 - Medical devices. Information to be supplied by the manufacturer, EN 60601-1 - 2006/A2:2021 - Medical Electrical Equipment - Part 1: General Requirements for Basic Safety and Essential Performance, EN 60601-1-2 - 2015/ A1:2021 - Medical Electrical Equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard:Electromagnetic Compatibility - Requirements and tests, EN 60601-1-6 - 2010+A1:2015+A2:2021 - Medical electrical equipment - Part 1-6: General requirements for basic safety and essential Performance - Collateral standard: Usability, EN 60601-1-11 - 2015+AD1:2020 - Medical Electrical Equipment -Part 1-11: General requirements for basic safety and essential performance -Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment, EN 62366 - 2015/A1:2020 - Medical devices - Application of usability engineering to medical devices, MEDDEV 2.7.1 - 2016 - GUIDELINES ON MEDICAL DEVICES - CLINICAL EVALUATION-A GUIDE FOR MANUFACTURERS AND NOTIFIED, ISO 27427 - 2023 - Anaesthetic and respiratory equipment - Nebulizing systems and components, EN ISO 10993-1 - 2020 - Biological evaluation of medical devices -Part 1: Evaluation and testing within a risk management process, EN ISO 10993-5 - 2009 - Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity, EN ISO 10993-10 - 2023 - Biological evaluation of medical devices - Part 10: Test for irritation and skin sensitization, EN ISO 10993-23 - 2021 - Biological evaluation of medical devices -Part 23: Tests for irritation, ISO 18562-2 - 2017 - Biocompatibility evaluation of breathing gas pathways in healthcare applications -Part 2: Tests for emissions of particulate matter, ISO 18562-3 - 2017 - Biocompatibility evaluation of breathing gas pathways in healthcare applications -Part 3: Tests for emissions of volatile organic compounds (VOCs), ISO 18562-4 - 2017 - Biocompatibility evaluation of breathing gas pathways in healthcare applications -Part 4: Tests for leachables in condensate.

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NORMES

Le produit est conforme aux normes suivantes : 93/42/EEC - 1993+Amd:2007 - Medical devices directives, EU 2017/745, EN ISO 14971 - 2019 - Medical devices-Application of risk management to Medical devices, EN ISO 13485 - 2016+A1:2021 - Medical devices - Quality management systems - Requirements for regulatory purposes, ISO 15223-1 - 2021 - Medical devices - Symbols to be used with information to be supplied by the manufacturer Part1:General Requirement, EN ISO 20417 - 2021 - Medical devices. Information to be supplied by the manufacturer, EN 60601-1 - 2006/A2:2021 - Medical Electrical Equipment - Part 1: General Requirements for Basic Safety and Essential Performance, EN 60601-1-2 - 2015/A1:2021 - Medical Electrical Equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard:Electromagnetic Compatibility - Requirements and tests, EN 60601-1-6 - 2010+A1:2015+A2:2021 - Medical electrical equipment - Part 1-6: General requirements for basic safety and essential Performance - Collateral standard: Usability, EN 60601-1-11 - 2015+AD1:2020 - Medical Electrical Equipment - Part 1-11:General requirements for basic safety and essential performance - Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment, EN 62366 - 2015/A1:2020 - Medical devices - Application of usability engineering to medical devices, MEDDEV 2.7.1 - 2016 - GUIDELINES ON MEDICAL DEVICES - CLINICAL EVALUATION-A GUIDE FOR MANUFACTURERS AND NOTIFIED, ISO 27427 - 2023 - Anaesthetic and respiratory equipment - Nebulizing systems and components, EN ISO 10993-1 - 2020 - Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process, EN ISO 10993-5 - 2009 - Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity, EN ISO 10993-10 - 2023 - Biological evaluation of medical devices - Part 10: Test for irritation and skin sensitization, EN ISO 10993-23 - 2021 - Biological evaluation of medical devices - Part 23: Tests for irritation, ISO 18562-2 - 2017 - Biocompatibility evaluation of breathing gas pathways in healthcare applications - Part 2: Tests for emissions of particulate matter, ISO 18562-3 - 2017 - Biocompatibility evaluation of breathing gas pathways in healthcare applications - Part 3: Tests for emissions of volatile organic compounds (VOCs), ISO 18562-4 - 2017 - Biocompatibility evaluation of breathing gas pathways in healthcare applications - Part 4: Tests for leachables in condensate.

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Το προϊόν ανταποκρίνεται στα παρακάτω πρότυπα: 93/42/EEC - 1993+Amd:2007 - Medical devices directives, EU 2017/745, EN ISO 14971 - 2019 - Medical devices-Application of risk management to Medical devices, EN ISO 13485 - 2016+A1:2021 - Medical devices - Quality management systems - Requirements for regulatory purposes, ISO 15223-1 - 2021 - Medical devices - Symbols to be used with information to be supplied by the manufacturer Part1:General Requirement, EN ISO 20417 - 2021 - Medical devices. Information to be supplied by the manufacturer, EN 60601-1 - 2006/A2:2021 - Medical Electrical Equipment - Part 1: General Requirements for Basic Safety and Essential Performance, EN 60601-1-2 - 2015/A1:2021 - Medical Electrical Equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic Compatibility - Requirements and tests, EN 60601-1-6 - 2010+A1:2015+A2:2021 - Medical electrical equipment - Part 1-6: General requirements for basic safety and essential Performance - Collateral standard: Usability, EN 60601-1-11 - 2015+AD1:2020 - Medical Electrical Equipment - Part 1-11: General requirements for basic safety and essential performance - Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment, EN 2015/A1:2020 - Medical devices - Application of usability engineering to medical devices, MEDDEV 2.7.1 - 2016 - GUIDELINES ON MEDICAL DEVICES - CLINICAL EVALUATION-A GUIDE FOR MANUFACTURERS AND NOTIFIED, ISO 27427 - 2023 - Anaesthetic and respiratory equipment - Nebulizing systems and components, EN ISO 10993-1 - 2020 - Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process, EN ISO 10993-5 - 2009 - Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity, EN ISO 10993-10 - 2023 - Biological evaluation of medical devices - Part 10: Test for irritation and skin sensitization, EN ISO 10993-23 - 2021 - Biological evaluation of medical devices - Part 23: Tests for irritation, ISO 18562-2 - 2017 - Biocompatibility evaluation of breathing gas pathways in healthcare applications - Part 2: Tests for emissions of particulate matter, ISO 18562-3 - 2017 - Biocompatibility evaluation of breathing gas pathways in healthcare applications - Part 3: Tests for emissions of volatile organic compounds (VOCs), ISO 18562-4 - 2017 - Biocompatibility evaluation of breathing gas pathways in healthcare applications - Part 4: Tests for leachables in condensate.

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STANDARDE

Produsul corespunde standardelor de mai jos: 93/42/EEC - 1993+Amd:2007 - Medical devices directives, EU 2017/745, EN ISO 14971 - 2019 - Medical devices-Application of risk management to Medical devices, EN ISO 13485 - 2016+A1:2021 - Medical devices - Quality management systems - Requirements for regulatory purposes, ISO 15223-1 - 2021 - Medical devices - Symbols to be used with information to be supplied by the manufacture Part1:General Requirement, EN ISO 20417 - 2021 - Medical devices. Information to be supplied by the manufacturer, EN 60601-1 - 2006/A2:2021 - Medical Electrical Equipment - Part 1: General Requirements for Basic Safety and Essential Performance, EN 60601-1-2 - 2015/A1:2021 - Medical Electrical Equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard:Electromagnetic Compatibility - Requirements and tests, EN 60601-1-6 - 2010+A1:2015+A2:2021 - Medical electrical equipment - Part 1-6: General requirements for basic safety and essential Performance - Collateral standard: Usability, EN 60601-1-11 - 2015+A1:2020 - Medical Electrical Equipment - Part 1-11:General requirements for basic safety and essential performance - Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment, EN 62366 - 2015/A1:2020 - Medical devices - Application of usability engineering to medical devices, MEDDEV 2.7.1 - 2016 - GUIDELINES ON MEDICAL DEVICES - CLINICAL EVALUATION-A GUIDE FOR MANUFACTURERS AND NOTIFIED, ISO 27427 - 2023 - Anaesthetic and respiratory equipment - Nebulizing systems and components, EN ISO 10993-1 - 2020 - Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process, EN ISO 10993-5 - 2009 - Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity, EN ISO 10993-10 - 2023 - Biological evaluation of medical devices - Part 10: Test for irritation and skin sensitization, EN ISO 10993-23 - 2021 - Biological evaluation of medical devices - Part 23: Tests for irritation, ISO 18562-2 - 2017 - Biocompatibility evaluation of breathing gas pathways in healthcare applications - Part 2: Tests for emissions of particulate matter, ISO 18562-3 - 2017 - Biocompatibility evaluation of breathing gas pathways in healthcare applications - Part 3: Tests for emissions of volatile organic compounds (VOCs), ISO 18562-4 - 2017 - Biocompatibility evaluation of breathing gas pathways in healthcare applications - Part 4: Tests for leachables in condensate.

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STANDARDS

Výrobek odpovídá níže uvedeným normám: 93/42/EEC - 1993+Amd:2007 - Medical devices directives, EU 2017/745, EN ISO 14971 - 2019 - Medical devices-Application of risk management to Medical devices, EN ISO 13485 - 2016+A1:2021 - Medical devices - Quality management systems - Requirements for regulatory purposes, ISO 15223-1 - 2021 - Medical devices – Symbols to be used with information to be supplied by the manufacturer Part1:General Requirement, EN ISO 20417 - 2021 - Medical devices. Information to be supplied by the manufacturer, EN 60601-1 - 2006/A2:2021 - Medical Electrical Equipment - Part 1: General Requirements for Basic Safety and Essential Performance, EN 60601-1-2 - 2015/A1:2021 - Medical Electrical Equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard:Electromagnetic Compatibility - Requirements and tests, EN 60601-1-6 - 2010+A1:2015+A2:2021 - Medical electrical equipment - Part 1-6: General requirements for basic safety and essential Performance - Collateral standard: Usability, EN 60601-1-11 - 2015+AD1:2020 - Medical Electrical Equipment - Part 1-11: General requirements for basic safety and essential performance - Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment, EN 62366 - 2015/A1:2020 - Medical devices - Application of usability engineering to medical devices, MEDDEV 2.7.1 - 2016 - GUIDELINES ON MEDICAL DEVICES - CLINICAL EVALUATION-A GUIDE FOR MANUFACTURERS AND NOTIFIED, ISO 27427 - 2023 - Anaesthetic and respiratory equipment - Nebulizing systems and components, EN ISO 10993-1 - 2020 - Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process, EN ISO 10993-5 - 2009 - Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity, EN ISO 10993-10 - 2023 - Biological evaluation of medical devices - Part 10: Test for irritation and skin sensitization, EN ISO 10993-23 - 2021 - Biological evaluation of medical devices - Part 23: Tests for irritation, ISO 18562-2 - 2017 - Biocompatibility evaluation of breathing gas pathways in healthcare applications – Part 2: Tests for emissions of particulate matter, ISO 18562-3 - 2017 - Biocompatibility evaluation of breathing gas pathways in healthcare applications – Part 3: Tests for emissions of volatile organic compounds (VOCs), ISO 18562-4 - 2017 - Biocompatibility evaluation of breathing gas pathways in healthcare applications – Part 4: Tests for leachables in condensate.

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Vyrobil: Vapo Healthcare Co., Ltd
Southern Unit of Third Floor, Building B, No.99 Yudai West Rd,
High Tech District, Kunshan, Suzhou, 215301, Jiangsu, P.R. China

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ŠTANDARDY

Výrobok zodpovedá nižšie uvedeným normám: 93/42/EEC - 1993+Amd:2007 - Medical devices directives, EU 2017/745, EN ISO 14971 - 2019 - Medical devices-Application of risk management to Medical devices, EN ISO 13485 - 2016+A1:2021 - Medical devices - Quality management systems - Requirements for regulatory purposes, ISO 15223-1 - 2021 - Medical devices – Symbols to be used with information to be supplied by the manufacturer Part1:General Requirement, EN ISO 20417 - 2021 - Medical devices. Information to be supplied by the manufacturer, EN 60601-1 - 2006/A2:2021 - Medical Electrical Equipment - Part 1: General Requirements for Basic Safety and Essential Performance, EN 60601-1-2 - 2015/A1:2021 - Medical Electrical Equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard:Electromagnetic Compatibility - Requirements and tests, EN 60601-1-6 - 2010+A1:2015+A2:2021 - Medical electrical equipment - Part 1-6: General requirements for basic safety and essential Performance - Collateral standard: Usability, EN 60601-1-11 - 2015+AD1:2020 - Medical Electrical Equipment - Part 1-11: General requirements for basic safety and essential performance - Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment, EN 62366 - 2015/A1:2020 - Medical devices - Application of usability engineering to medical devices, MEDDEV 2.7.1 - 2016 - GUIDELINES ON MEDICAL DEVICES - CLINICAL EVALUATION-A GUIDE FOR MANUFACTURERS AND NOTIFIED, ISO 27427 - 2023 - Anaesthetic and respiratory equipment - Nebulizing systems and components, EN ISO 10993-1 - 2020 - Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process, EN ISO 10993-5 - 2009 - Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity, EN ISO 10993-10 - 2023 - Biological evaluation of medical devices - Part 10: Test for irritation and skin sensitization, EN ISO 10993-23 - 2021 - Biological evaluation of medical devices - Part 23: Tests for irritation, ISO 18562-2 - 2017 - Biocompatibility evaluation of breathing gas pathways in healthcare applications – Part 2: Tests for emissions of particulate matter, ISO 18562-3 - 2017 - Biocompatibility evaluation of breathing gas pathways in healthcare applications – Part 3: Tests for emissions of volatile organic compounds (VOCs), ISO 18562-4 - 2017 - Biocompatibility evaluation of breathing gas pathways in healthcare applications – Part 4: Tests for leachables in condensate.

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STANDARDS

A termék megfelel az alábbi szabványoknak: 93/42/EEC - 1993+Amd:2007 - Medical devices directives, EU 2017/745, EN ISO 14971 - 2019 - Medical devices-Application of risk management to Medical devices, EN ISO 13485 - 2016+A1:2021 - Medical devices - Quality management systems - Requirements for regulatory purposes, ISO 15223-1 - 2021 - Medical devices - Symbols to be used with information to be supplied by the manufacturer Part1:General Requirement, EN ISO 20417 - 2021 - Medical devices. Information to be supplied by the manufacturer, EN 60601-1 - 2006/A2:2021 - Medical Electrical Equipment - Part 1: General Requirements for Basic Safety and Essential Performance, EN 60601-1-2 - 2015/ A1:2021 - Medical Electrical Equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard:Electromagnetic Compatibility - Requirements and tests, EN 60601-1-6 - 2010+A1:2015+A2:2021 - Medical electrical equipment - Part 1-6: General requirements for basic safety and essential Performance - Collateral standard: Usability, EN 60601-1-11 - 2015+AD1:2020 - Medical Electrical Equipment - Part 1-11: General requirements for basic safety and essential performance - Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment, EN 62366 - 2015/A1:2020 - Medical devices - Application of usability engineering to medical devices, MEDDEV 2.7.1 - 2016 - GUIDELINES ON MEDICAL DEVICES - CLINICAL EVALUATION-A GUIDE FOR MANUFACTURERS AND NOTIFIED, ISO 27427 - 2023 - Anaesthetic and respiratory equipment - Nebulizing systems and components, EN ISO 10993-1 - 2020 - Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process, EN ISO 10993-5 - 2009 - Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity, EN ISO 10993-10 - 2023 - Biological evaluation of medical devices - Part 10: Test for irritation and skin sensitization, EN ISO 10993-23 - 2021 - Biological evaluation of medical devices - Part 23: Tests for irritation, ISO 18562-2 - 2017 - Biocompatibility evaluation of breathing gas pathways in healthcare applications - Part 2: Tests for emissions of particulate matter, ISO 18562-3 - 2017 - Biocompatibility evaluation of breathing gas pathways in healthcare applications - Part 3: Tests for emissions of volatile organic compounds (VOCs), ISO 18562-4 - 2017 - Biocompatibility evaluation of breathing gas pathways in healthcare applications - Part 4: Tests for leachables in condensate.

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STANDARTI

Izstrādājums atbilst turpmāk minētajiem standartiem: 93/42/EEC - 1993+Amd:2007 - Medical devices directives, EU 2017/745, EN ISO 14971 - 2019 - Medical devices-Application of risk management to Medical devices, EN ISO 13485 - 2016+A1:2021 - Medical devices - Quality management systems - Requirements for regulatory purposes, ISO 15223-1 - 2021 - Medical devices — Symbols to be used with information to be supplied by the manufacturer Part1:General Requirement, EN ISO 20417 - 2021 - Medical devices. Information to be supplied by the manufacturer, EN 60601-1 - 2006/A2:2021 - Medical Electrical Equipment - Part 1: General Requirements for Basic Safety and Essential Performance, EN 60601-1-2 - 2015/A1:2021 - Medical Electrical Equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic Compatibility - Requirements and tests, EN 60601-1-6 - 2010+A1:2015+A2:2021 - Medical electrical equipment- Part 1-6: General requirements for basic safety and essential Performance - Collateral standard: Usability, EN 60601-1-11 - 2015+AD1:2020 - Medical Electrical Equipment - Part 1-11: General requirements for basic safety and essential performance - Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment, EN 62366 - 2015/A1:2020 - Medical devices - Application of usability engineering to medical devices, MEDDEV 2.7.1 - 2016 - GUIDELINES ON MEDICAL DEVICES - CLINICAL EVALUATION-A GUIDE FOR MANUFACTURERS AND NOTIFIED, ISO 27427 - 2023 - Anaesthetic and respiratory equipment - Nebulizing systems and components, EN ISO 10993-1 - 2020 - Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process, EN ISO 10993-5 - 2009 - Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity, EN ISO 10993-10 - 2023 - Biological evaluation of medical devices - Part 10: Test for irritation and skin sensitization, EN ISO 10993-23 - 2021 - Biological evaluation of medical devices - Part 23: Tests for irritation, ISO 18562-2 - 2017 - Biocompatibility evaluation of breathing gas pathways in healthcare applications — Part 2: Tests for emissions of particulate matter, ISO 18562-3 - 2017 - Biocompatibility evaluation of breathing gas pathways in healthcare applications — Part 3: Tests for emissions of volatile organic compounds (VOCs), ISO 18562-4 - 2017 - Biocompatibility evaluation of breathing gas pathways in healthcare applications — Part 4: Tests for leachables in condensate.

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Продуктът отговаря на следните стандарти: 93/42/EEC - 1993+Amd:2007 - Medical devices directives, EU 2017/745, EN ISO 14971 - 2019 - Medical devices-Application of risk management to Medical devices, EN ISO 13485 - 2016+A1:2021 - Medical devices - Quality management systems - Requirements for regulatory purposes, ISO 15223-1 - 2021 - Medical devices – Symbols to be used with information to be supplied by the manufacturer Part1:General Requirement, EN ISO 20417 - 2021 - Medical devices. Information to be supplied by the manufacturer, EN 60601-1 - 2006/A2:2021 - Medical Electrical Equipment - Part 1: General Requirements for Basic Safety and Essential Performance, EN 60601-1-2 - 2015/A1:2021 - Medical Electrical Equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard:Electromagnetic Compatibility - Requirements and tests, EN 60601-1-6 - 2010+A1:2015+A2:2021 - Medical electrical equipment - Part 1-6: General requirements for basic safety and essential Performance - Collateral standard: Usability, EN 60601-1-11 - 2015+AD1:2020 - Medical Electrical Equipment - Part 1-11: General requirements for basic safety and essential performance - Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment, EN 62366 - 2015/A1:2020 - Medical devices - Application of usability engineering to medical devices, MEDDEV 2.7.1 - 2016 - GUIDELINES ON MEDICAL DEVICES - CLINICAL EVALUATION-A GUIDE FOR MANUFACTURERS AND NOTIFIED, ISO 27427 - 2023 - Anaesthetic and respiratory equipment - Nebulizing systems and components, EN ISO 10993-1 - 2020 - Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process, EN ISO 10993-5 - 2009 - Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity, EN ISO 10993-10 - 2023 - Biological evaluation of medical devices - Part 10: Test for irritation and skin sensitization, EN ISO 10993-23 - 2021 - Biological evaluation of medical devices - Part 23: Tests for irritation, ISO 18562-2 - 2017 - Biocompatibility evaluation of breathing gas pathways in healthcare applications – Part 2: Tests for emissions of particulate matter, ISO 18562-3 - 2017 - Biocompatibility evaluation of breathing gas pathways in healthcare applications – Part 3: Tests for emissions of volatile organic compounds (VOCs), ISO 18562-4 - 2017 - Biocompatibility evaluation of breathing gas pathways in healthcare applications – Part 4: Tests for leachables in condensate.

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STANDARTAI

Gaminys atitinka toliau nurodytus standartus: 93/42/EEC - 1993+Amd:2007 - Medical devices directives, EU 2017/745, EN ISO 14971 - 2019 - Medical devices-Application of risk management to Medical devices, EN ISO 13485 - 2016+A1:2021 - Medical devices - Quality management systems - Requirements for regulatory purposes, ISO 15223-1 - 2021 - Medical devices – Symbols to be used with information to be supplied by the manufacturer Part1:General Requirement, EN ISO 20417 - 2021 - Medical devices. Information to be supplied by the manufacturer, EN 60601-1 - 2006/A2:2021 - Medical Electrical Equipment - Part 1: General Requirements for Basic Safety and Essential Performance, EN 60601-1-2 - 2015/A1:2021 - Medical Electrical Equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard:Electromagnetic Compatibility - Requirements and tests, EN 60601-1-6 - 2010+A1:2015+A2:2021 - Medical electrical equipment - Part 1-6: General requirements for basic safety and essential Performance - Collateral standard: Usability, EN 60601-1-11 - 2015+AD1:2020 - Medical Electrical Equipment - Part 1-11: General requirements for basic safety and essential performance - Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment, EN 62366 - 2015/A1:2020 - Medical devices - Application of usability engineering to medical devices, MEDDEV 2.7.1 - 2016 - GUIDELINES ON MEDICAL DEVICES - CLINICAL EVALUATION-A GUIDE FOR MANUFACTURERS AND NOTIFIED, ISO 27427 - 2023 - Anaesthetic and respiratory equipment - Nebulizing systems and components, EN ISO 10993-1 - 2020 - Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process, EN ISO 10993-5 - 2009 - Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity, EN ISO 10993-10 - 2023 - Biological evaluation of medical devices - Part 10: Test for irritation and skin sensitization, EN ISO 10993-23 - 2021 - Biological evaluation of medical devices - Part 23: Tests for irritation, ISO 18562-2 - 2017 - Biocompatibility evaluation of breathing gas pathways in healthcare applications – Part 2: Tests for emissions of particulate matter, ISO 18562-3 - 2017 - Biocompatibility evaluation of breathing gas pathways in healthcare applications – Part 3: Tests for emissions of volatile organic compounds (VOCs), ISO 18562-4 - 2017 - Biocompatibility evaluation of breathing gas pathways in healthcare applications – Part 4: Tests for leachables in condensate.

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Parengė: Vapo Healthcare Co., Ltd
Southern Unit of Third Floor, Building B, No.99 Yudai West Rd,
High Tech District, Kunshan, Suzhou, 215301, Jiangsu, P.R. China

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STANDARDID

Toode vastab allpool esitatud standarditele: 93/42/EEC - 1993+Amd:2007 - Medical devices directives, EU 2017/745, EN ISO 14971 - 2019 - Medical devices-Application of risk management to Medical devices, EN ISO 13485 - 2016+A1:2021 - Medical devices - Quality management systems - Requirements for regulatory purposes, ISO 15223-1 - 2021 - Medical devices – Symbols to be used with information to be supplied by the manufacturer Part 1: General Requirement, EN ISO 20417 - 2021 - Medical devices. Information to be supplied by the manufacturer, EN 60601-1 - 2006/A2:2021 - Medical Electrical Equipment - Part 1: General Requirements for Basic Safety and Essential Performance, EN 60601-1-2 - 2015/A1:2021 - Medical Electrical Equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic Compatibility - Requirements and tests, EN 60601-1-6 - 2010+A1:2015+A2:2021 - Medical electrical equipment - Part 1-6: General requirements for basic safety and essential Performance - Collateral standard: Usability, EN 60601-1-11 - 2015+AD1:2020 - Medical Electrical Equipment - Part 1-11: General requirements for basic safety and essential performance - Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment, EN 62366 - 2015/A1:2020 - Medical devices - Application of usability engineering to medical devices, MEDDEV 2.7.1 - 2016 - GUIDELINES ON MEDICAL DEVICES - CLINICAL EVALUATION-A GUIDE FOR MANUFACTURERS AND NOTIFIED, ISO 27427 - 2023 - Anaesthetic and respiratory equipment - Nebulizing systems and components, EN ISO 10993-1 - 2020 - Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process, EN ISO 10993-5 - 2009 - Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity, EN ISO 10993-10 - 2023 - Biological evaluation of medical devices - Part 10: Test for irritation and skin sensitization, EN ISO 10993-23 - 2021 - Biological evaluation of medical devices - Part 23: Tests for irritation, ISO 18562-2 - 2017 - Biocompatibility evaluation of breathing gas pathways in healthcare applications – Part 2: Tests for emissions of particulate matter, ISO 18562-3 - 2017 - Biocompatibility evaluation of breathing gas pathways in healthcare applications – Part 3: Tests for emissions of volatile organic compounds (VOCs), ISO 18562-4 - 2017 - Biocompatibility evaluation of breathing gas pathways in healthcare applications – Part 4: Tests for leachables in condensate.

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Tootja: Vapo Healthcare Co., Ltd
Jagab: Laica S.p.A.
LAICA S.p.A.