IT Italiano

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ELECTRONIC PERSONAL SCALE

Instructions and warranty

Your new bathroom scale is extremely reliable and precise, it does not require any sort of maintenance. To improve the reading, weigh yourself twice, and if the weights are different, your weight is the middle of the two readings. Always weigh yourself on the same scale, at the same time of day, preferably undressed and before breakfast. Caution! The surface of the scale can be slippery, if wet.

IT IS IMPORTANT TO READ THE INSTRUCTIONS AND WARNINGS IN THIS HANDBOOK BEFORE USE AND KEEP IT IN A SAFE PLACE.

THE BODY WEIGHT

Checking the weight means performing an **energetic balance**: such balance results from the ratio between the quantity of ingested nutrients, i.e. the ingested calories, and the actually consumed calories. The daily necessary quantity of energy depends on the metabolism, on the life-style and on the type of activities performed by each single person.

Controlling one's own weight is essential in order to protect health effectively: eating well means being healthy. It has been proved that most of diseases, many of which extremely dangerous, are associated to unbalanced eating habits not only in terms of quality but also in terms of quantity. Among the most known, the following ones should be mentioned:

- ischemic heart diseases;
- cerebrovascular diseases;
- diabetes mellitus;
- intestine tumors;
- etc. ...

All these diseases are present above all among obese persons: therefore, controlling the weight together with controlling the quality of the nutrients becomes one of the essential factors to prevent diseases. Healthy eating means selecting the foodstuffs that supply the body with all the necessary nutriment to grow and keep the physiological functions and the biochemical processes balanced.

Hints for a healthy feeding and to control weight:

- To follow a balanced diet it is necessary not to exceed in the consumption of the same types of food neglecting the essential elements contained in other types of food, as well as not eating more than needed.
- A healthy diet includes suitable quantities of proteins, carbohydrates, fats, fibers, vitamins and minerals, as well as, obviously, water.
- Having a complete meal only once a day is a bad habit. It is better to have many small meals instead of just an abundant one.
- Three meals a day and a snack in the middle of the morning or of the afternoon represent the basis for a correct feeding.
- It is suggested to sit down and relax when eating. Select a comfortable environment can be relevant.
- Do not eat when you are focused on doing other activities, i.e. when you are watching TV or working, in order to avoid eating too much without even noticing it.
- To support the digestion process it is recommended to eat slowly and chew well, enjoying the color, taste and consistency of the food.
- Combining a correct feeding to physical activity is an essential condition to allow increasing the metabolism, consuming more calories and remaining in good shape.
- It is important to remember that with age human body seems to need less energy due to still not so clear reasons (maybe due to the fact that one leads a less active life); however, many aged people do not reduce their consumption of food and tend to grow fat.
- Getting used to follow a balanced diet is surely of help in increasing the quality of life, promoting health and preventing diseases, as well as in increasing vitality and remaining in good shape!

Caution!!!!

- Never follow diets autonomously, always refer to a doctor or dietician.**
- Self-measurement means control, not diagnosis or treatment. Unusual values have always to be discussed with one's own physician.**
- Under no circumstances must you change the dose of medicines prescribed by your Doctor.**

ISTRUCTIONS FOR USE

- Place the scale on a flat, hard surface. If the floor is loose or uneven (e.g. mats, carpets, linoleum) the exactness of your weight is not guaranteed.**
- Place the scale in a well-lit area. When  flashes on the display, the scale is ready to use.
- Select the weighing unit of measure (kg, lb or st) with the key located on the bottom of the scale. The unit of measurement can be selected when  flashes on the display.
- Step onto the scale with your feet parallel and your weight equally distributed (see illustrations). Warning: the scale switches on once step on. Stand still while the scale measures your weight. Do not lean on anything while being weighed as this could alter the reading.**
- Step off the scale; the display screen will continue to show your weight for a few seconds then the scale will turn itself off automatically.

WARNING:

- In case of overweight, the display screen shows "Err", remove the weight from the scale to avoid any damage of the mechanism.
- When on the display there is no  it is necessary to expose the scale to a sufficiently lightened environment for at least 1 minute.

MAINTENANCE

- Handle your scale with care as it is a precision instrument. Therefore, avoid to expose it to extreme temperatures, damp, knocks, dust and direct sunlight.
- Do not drop, seat, cause or knock on it.
- Keep the scale away from sources of heat.
- Do not jump on the scale.**
- Clean the platform with a damp cloth.**
- Clean the bottom of the scale with just a soft dry cloth.**
- To clean the scale never use chemical or abrasive products.**
- Take care that no liquids penetrate into the scale, as it is not waterproof.
- Do not leave or use the scale in very damp surroundings, humidity above 85% or contact with water and other liquid parts could cause corrosion and jeopardize the precision of the scale.

TECHNICAL DATAS

LCD display
Maximum tolerance: +/-1% +1 division
Working conditions: +10°C +40°C; RH 85%
Storage conditions: +10°C +50°C; RH 85%

i¹ Important! Read the instructions carefully

This product is meant for domestic use. The conformity, certified by the CE-marking on the device, refers to the EEC Directive 2004/108 on electromagnetic compatibility. This product is an electronic device that was tested to assure, according to the current technical knowledge, that it does not interfere with other nearby located equipment (electromagnetic compatibility) and that it is safe if used according to

GB English

ELECTRONIC PERSONAL SCALE

Instructions and warranty

the advices given in the instructions for use. In case of operation anomalies, do not keep using the device and, if necessary, refer to the manufacturer directly.

DISPOSAL PROCEDURE (Dir. 2012/19/Eu-WEEE)

The symbol on the bottom of the device indicates the separated collection of electric and electronic equipment.

At the end of life of the device, do not dispose it as mixed solid municipal waste, but dispose it referring to a specific collection centre located in your area or returning it to the distributor, when buying a new device of the same type to be used with the same functions.

If the appliance to be disposed of is less than 25 cm, it can be returned to a retail location that is over 400 m² without having to purchase a new, similar device. This procedure of separated collection of electric and electronic devices is carried out forecasting a European environmental policy aiming at safeguarding, protecting and improving environment quality, as well as avoiding potential effects on human health due to the presence of hazardous substances in such equipment or to an improper use of the same or of parts of the same. Caution! The wrong disposal of electric and electronic equipment may involve sanctions.

WARRANTY

This equipment is covered by a 2-year warranty starting from the purchase date, which has to be proven by means of **the stamp or signature of the dealer and by the fiscal payment receipt that shall be kept enclosed to this document**. This warranty period complies with the laws in force and applies only in case the consumer is an individual. Laica products are designed for domestic use only; therefore, their use for commercial purposes is forbidden. The warranty covers only production defects and does not apply, if the damage is caused by accidental events, misuse, negligence or improper use of the product. Use the supplied accessories only; the use of other accessories invalidate the warranty cover.

Do not open the equipment for any reason whatsoever; the opening or the tampering of the equipment, definitely invalidate the warranty. Warranty does not apply to parts subject to tear and wear along the use, as well as to batteries, when they are included.

Warranty decays after 2 years from the purchase date. In this case, technical service interventions will be carried out against payment. Information on technical service interventions, under warranty terms or against payment, may be requested by writing at info@laica.com. Repair or replacement of products covered by the warranty terms are free of charge. In case of faults, refer to the dealer; DO NOT send anything directly to LAICA.

All interventions covered by warranty (including product replacement or the replacement of any of its parts) shall not be meant to extend the duration of the original warranty period of the replaced product.

The manufacturer cannot be made liable for any damages that may be caused, directly or indirectly, to people, things and pets resulting from the inobservance of the provisions contained in the suitable user manual, especially with reference to the instructions related to installation, use and maintenance of the equipment.

Being constantly involved in improving its own products, Laica reserves the right to modify its product, totally or partially, with reference to any production need, without that this involves any responsibility whatsoever by Laica or its dealers.

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