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NORME

Il dispositivo è conforme alle norme: ISO80601-2-56 Medical electrical equipment - Part 2-56:Particular requirements for basic safety and essential performance of clinical thermometers for body temperature measurement, EN 60601-1-11 Medical electrical equipment - Part 1-11: General requirements for basic safety and essential performance - Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment and complies with the requirements of EN 60601-1-2(EMC), IEC/EN60601-1(Safety) standards.

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NON-CONTACT FOREHEAD INFRARED THERMOMETER INSTRUCTIONS AND WARRANTY

Dear Customer, Laica would like to thank you for choosing this product, designed according to criteria of reliability and quality for complete satisfaction.

IMPORTANT READ CAREFULLY BEFORE USING - KEEP IN A SAFE PLACE FOR FUTURE REFERENCE

The instruction manual must be considered an integral part of the product and must be kept for its entire life. If the product is transferred to another owner, its documentation must also be transferred in its entirety.

To ensure the safe and correct use of the product, the user must carefully read the instructions and warnings contained in the manual insofar as they provide important information regarding safety, user and maintenance instructions.

Should the instruction manual get mislaid or you require additional information or clarifications, please fill in the relevant form on the website: <https://www.laica.it/> in the FAQ and Support section.

The LAICA infrared thermometer measures the heat emitted from the forehead over the temporal artery to indicate the patient's body temperature. The innovative remote mode for front detection makes it ideal for measurements on children, even during sleep. It is very easy to use and its purpose is to detect body temperature and the temperature of liquids. With large LCD display. N.B.: the body temperature shown on the display is equivalent to the oral temperature.

KEY TO SYMBOLS

Warning

Prohibition

Attention! Read the user instructions carefully

Symbol of "type BF applied parts"

C €0425

In compliance with European legislation on medical devices

Manufacturer

Production batch number

IP22: Degree of protection of enclosures for electrical equipment, where the first digit indicates the degree of protection against the penetration of solid foreign bodies (from 0 to 6) and the second digit the degree of protection against the penetration of liquids (from 0 to 8).

SAFETY WARNINGS

- Before using the product, check that the appliance is intact without visible damage. If in doubt, do not use the product and contact your dealer.
- This device must be used only for its intended purpose and according to the user instructions. All other types of use shall be considered inappropriate and therefore hazardous. The manufacturer shall not be deemed liable for any damage caused by inappropriate or incorrect use.
- People with physical, sensory or mental disabilities or unskilled people, may use and perform maintenance on the product only under adult supervision. Children must not play with the device.

- This appliance can be used autonomously by children of 12 years of age or older (able to read and understand the instructions for use).
- Treat the product with care, protect it from impacts, extreme variations in temperature, humidity, dust, direct sunlight and sources of heat.
- In the event of a failure and/or malfunction, turn off the device without tampering with it. Always contact the retailer for repairs.
- Ensure hands are dry when using the on and off key. NEVER submerge the product in water or any other liquids.
- Device for household use, which must not be used to measure temperature outdoors.

ATTENTION! BEFORE USING THIS APPLIANCE

- 1) **Self-measurement means checking, not diagnosis or treatment. Unusual values must always be discussed with your doctor. Under no circumstances should you alter the dosages of any medication prescribed by your doctor.**
- 2) This thermometer is a sensitive instrument. Handle with care and do not expose it to mechanical shocks. Do not squeeze, bend, drop or take apart the thermometer.
- 3) Do not touch with your fingers and do not blow on the probe; do not disassemble the thermometer.
- 4) Do not use this thermometer to detect rectal, oral or underarm temperature. It must be used to measure the forehead temperature at a distance of 1-3 cm.
- 5) Measurement accuracy can be compromised by: pronounced sweating of the forehead, taking vasoconstricting medication and skin irritations. Holding the thermometer in your hand for too long can cause a false recording of body temperature.
- 6) Do not measure the temperature on wounds, abrasions or scars.

- 7) Remove any hair from the forehead at least 10 minutes before taking the measurement. Clean the forehead where the probe will be pointed before each recording.
- 8) **Measure the temperature far from direct sunlight, wind, air conditioning and sources of heat such as stoves and fireplaces.**
- 9) Do not measure the temperature immediately after breastfeeding a newborn or for thirty minutes after physical exercise or swimming.
- 10) The appliance may provide inaccurate readings when used in temperature or humidity conditions outside the limits indicated in the "Technical Specifications" paragraph.
- 11) Do not use near strong magnetic fields and keep away from radio systems or mobile phones (for more information on interferences see paragraph "Electromagnetic compatibility").
- 12) The appliance is not suitable for use in the presence of inflammable anaesthetic mixes with air, oxygen, or nitrous oxide.

SAFE USE OF BATTERIES

- Remove the batteries if the device is not used for long periods of time and keep them in a well-ventilated and dry place, at room temperature.
- Do NOT use rechargeable batteries.
- NEVER expose the batteries to sources of heat or direct sunlight. Failure to observe these instructions may result in damage to and/or the explosion of the batteries.
- Do NOT throw the batteries into an open flame.
- Batteries must always be removed or replaced by adults only.
- Keep batteries out of reach of children: the ingestion of batteries can cause serious harm and death. If a battery is swallowed, seek urgent medical care.
- The acid in batteries is corrosive. Avoid contact with the skin, eyes or garments.

FEVER

Fever should be considered a symptom and not a disease. In general, it is a sign that the body is fighting an infection. In fact, it is a defence mechanism of the body characterised by an increase in body temperature above **37°C** and increased heart and breathing rates. Normal temperature within the body is about 37.5°C and this is in fact the value measured with the rectal method. The temperature measured in the mouth is about 0.5°C less (37°C) while the underarm recording is 1°C less (36.5°C). From this we can deduce the importance of a correct recording of the temperature that should be carried out by valid means, and if possible not at a time when the body temperature may be increased for physiological reasons (for example after meals).

In a healthy person the temperature is affected by many factors:

- individual character of the person (individual metabolism);
- age (infants and toddlers' body temperature is higher and decreases with age). In children, the temperature varies with greater intensity, speed and frequency;
- clothing;
- outside air temperature;
- time of day (in the morning body temperature is lower and increases during the day);
- motor and mental activities carried out;
- measurement method;
- phase of the menstrual cycle.

It should be noted that the thermometer must already be in the room in which you will perform the measurement, about half an hour before taking the measurement, while the person must be in the room at least 10 minutes beforehand. The table below shows a list of "normal" average temperatures: **it is however advisable to become familiar with recognising your normal temperature when you are physically well in order to be able to interpret the different values recorded.**

METHOD OF MEASURING AVERAGE

UNDERARM

ORAL

RECTAL

EAR

FOREHEAD

NORMAL TEMPERATURE

36,5°C

37°C

37,5°C

37°C

37°C

Each type of thermometer lends itself to the measurement of body temperature at a specific point of the body: the forehead thermometer only for forehead measurement; the ear thermometer only for measurement in the ear; the digital thermometer for traditional rectal, axillary and oral measurement. Depending on the point where the measurement is made, a different temperature value is obtained. The variation can be between 0.2-1°C. You can not compare the temperatures measured with different measuring methods. It is therefore necessary to consider in case of self-diagnosis or to report to your doctor what type of thermometer and at which point of the body the temperature has been measured.

PRODUCT DESCRIPTION (see fig.1)

- 1) Measurement probe
- 2) "I" key
- 3) "SET" key
- 4) "O/I" key
- 5) Battery compartment
- 6) Batteries
- 7) LCD DISPLAY
 - a. Indication of forehead temperature measurement
 - b. Indication of liquid temperature measurement
 - c. Temperature recorded
 - d. Indication of the measurement result

- e. Unit of measure °C/°F
- f. Memory number
- g. Low battery indicator
- h. Date and time
- i. AM and PM indicator

TECHNICAL CHARACTERISTICS

- Product name: infra-red thermometer
- Trade name: TH1000
- Body temperature measurement range: from 34°C to 43°C (93.2°F - 109.4°F)
- Liquid temperature measurement range: from 0°C to 100°C (32°F - 212°F)
- Measurement zone: forehead
- Body reference zone: mouth (the thermometer converts the temperature recorded on the forehead and displays the equivalent oral temperature)
- Display resolution: 0.1°C/0.1°F
- Remote measurement: 1 – 3 cm
- Batteries: 2x1.5V alkaline AAA
- Switch off: automatic, after 60 seconds of non-use
- Battery autonomy: about 6000 readings (turning the thermometer on and off after each measurement) / 1 year
- Operating conditions: from 10°C to 40°C (50°F - 104°F), relative humidity 15%-85% non-condensing, atmospheric pressure: 700 – 1060 hPa
- Transport and storage conditions: from -25°C to 55°C (-13°F - 131°F), relative humidity 15%-85% non-condensing, atmospheric pressure: 700 – 1060 hPa
- Tolerance in lab for body temperature:
 - with ambient temperature between 15 and 35°C (59-95°F):

- from 35.5 to 42°C (95.9-107.6°F) : +/-0.2°C (0.4°F)
- with ambient temperature less than 15°C and higher than 35°C: +/-0.3°C (0.5°F)
- < 34.9°C (94.8°F) and > 42.1°C (107.7°F): +/- 0.3°C (0.5°F)
- Tolerance in lab for liquid temperature: +/-4% or +/-2°C(4°F)
- Clinical tolerance:
 - Clinical bias: -0.009°C (-0.16°F)
 - Clinical repeatability: 0.13°C (0.23°F)069
 - Acceptability limit: 0.87°C (1.57°F)

INSTRUCTIONS FOR USE

INSERTION/REPLACEMENT OF THE BATTERY

This thermometer operates with 2 x 1.5V AAA alkaline batteries. When first used, and when the display shows the battery symbol , proceed with insertion and/or battery replacement. Open the battery compartment by sliding the cover in the direction of the arrow, insert the batteries according to the pole markings and close the cover. Dispose of used batteries as indicated in the paragraph "Disposal procedure".

SETTINGS

• DATE AND TIME

- 1) When the device is switched off, press the "O/I" key for about 3 seconds: the thermometer will turn on and the 24 or 12 hour time format will show on the display. Use the "I" key to adjust the date and confirm using the "O/I" key.
- 2) Set the hours, minutes, year, month and day in order. Adjust the data using the "I" key and confirm using the "O/I" key. Hold the "I" key pressed to speed up the number scrolling. The display will show "OFF" and turn off. The date/time will need to be reset whenever the batteries are replaced.

• UNIT OF MEASUREMENT °C / °F

This thermometer can measure temperature in Celsius and Fahrenheit.

With the thermometer turned off, press the "0" key for about 3 seconds: the display will show the unit of measure (°C or °F). To change the unit of measure press the "0" key. Press the "O/I" key to turn the device off.

TEMPERATURE READING:

• FOREHEAD

Attention!

In order to ensure accurate recording of the temperature, clean the probe in its entirety with a cotton swab moistened with alcohol after each use. Make sure that the probe is completely dry before use.

AN INADEQUATE CLEANING OF THE PROBE CAN PROVIDE FALSE READINGS

- 1) Clean and dry the forehead where you will be positioning the probe.
- 2) Turn on the thermometer by pressing the "O/I" key.
The thermometer will emit an acoustic signal and all the function symbols will appear on the display (self-test). It will then emit two short acoustic signals: it is now possible to take the measurement.
- 3) Point the probe at the centre of the forehead just above the space between the eyebrows at a distance of approximately 1-3 cm and press the "0" key. The device will emit an acoustic signal and the display will show the recorded measurement, then emit two short acoustic signals. Do not move the thermometer before the two acoustic signals have sounded, confirming the recorded measurement. The display shows the "☺" symbol when the body temperature is less than 37.8°C (100°F) and the "☹" symbol when it is higher. If the temperature is higher than 37.8°C the thermometer emits a long acoustic signal followed by a series of short acoustic signals.

 The temperature shown on the display is the oral temperature.

- 4) The device automatically turns off after about 60 seconds, otherwise it can be turned off by pressing the "O/I" key.

Wait at least one second between one measurement and the next.

• SURFACES AND LIQUIDS

This thermometer can record the temperature of warm surfaces or liquids (such as milk, water,...).

1. Turn on the thermometer by pressing the "O/I" key. The thermometer will emit an acoustic signal and all the function symbols will appear on the display (self-test). It will then emit two short acoustic signals.
2. Hold the "SET" key pressed for about 3 seconds: the display shows the "🏠" symbol.
3. To measure the temperature, hold the thermometer facing the interested surface at a distance of about 1-3 cm and press the "0" key. The device will emit an acoustic signal and the display will show the recorded measurement, then emit two short acoustic signals.
Do not move the thermometer before the two acoustic signals have sounded, confirming the recorded measurement.
4. The device automatically turns off after about 60 seconds, otherwise it can be turned off by pressing the "O/I" key.

 **Never immerse the probe in any kind of liquid.**

MEMORY FUNCTION

This device saves the last 10 measurements recorded for each type of measurement (body and liquid) to evaluate any variations. Readings are automatically stored. Once 10 readings have been exceeded, the oldest data is automatically deleted.

- 1) Select the measuring mode to be consulted (body "👤" or surfaces and liquids "🏠") as described in the sections "Temperature reading – Forehead – Surfaces and liquids"

- 2) To call up the saved measurements, hold the "O/I" key pressed for about 3 seconds with the device on: an "M" will appear in the top right of the display with the last recorded measurement.
- 3) Press the "O/I" key again, keep pressing it to scroll through the stored data: number 1 refers to the most recent measurement, number 10 to the oldest measurement.
- 4) The device automatically turns off after about 60 seconds, otherwise it can be turned off by holding the "O/I" key pressed for about 3 seconds.

MAINTENANCE

Before and/or after each use, clean the probe of the thermometer with a cotton swab moistened with alcohol to ensure accurate readings. The probe is the most delicate part of the thermometer, pay particular attention to it during cleaning. Allow the probe to dry completely for at least 30 minutes.

The body of the thermometer should be cleaned with a soft, dry cloth. THE BODY OF THE THERMOMETER MUST NOT BE WASHED WITH WATER.

Never use abrasive products and do not immerse the thermometer in water or other liquids. Keep the thermometer in a clean, dry place, away from direct sunlight. It is advisable to check the performance of the appliance every 2 years or after a repair.

Contact the Laica help service (activity excluded from the warranty).

TROUBLESHOOTING

Problem	Possible cause	Solution
The device does not turn on.	The batteries are worn.	Replace the batteries.
	The batteries have not been installed correctly.	Check the polarity of the batteries during insertion.
The display shows 	The batteries are flat.	Replace the batteries.

Problem	Possible cause	Solution
The display shows "Hi" or "Lo".	The measured temperature does not fall within the thermometer measurement range: • body measurement: from 34°C to 43°C, • surface and liquid measurement: from 0°C to 100°C.	
	The temperature of a surface or a liquid is being measured without changing the thermometer measurement setting.	To measure the temperature of surfaces or liquids it is necessary to select the "🏠" measuring mode.
The display shows "Lo".	During forehead detection, the measurement distance is greater than 3 cm.	Approach the thermometer to the forehead during measurement at a distance of between 1 and 3 cm.
	The probe is dirty.	Clean the probe after each use as described in the "Maintenance" paragraph.
The display shows "Er1"	The measurement was made during the self-test.	To measure the temperature, wait for the display to show "👤" / "🏠"

Problem	Possible cause	Solution
The display shows "Er2"	The room temperature is not between 10°C and 40°C (50°F and 104°F).	Use the thermometer respecting the environmental conditions of operation (read the paragraph "Technical specifications").
The display shows "Er3"	During measurement the thermometer was positioned incorrectly.	Re-read the sections "Temperature reading - Forehead", "Temperature reading - Surfaces and liquids".
The display shows "Er4".	When measuring the temperature of surfaces or liquids, the thermometer is in an environment where the temperature rapidly fluctuates.	The thermometer must be positioned in the room where the measurement will be made at least 30 minutes beforehand, with a temperature between 10°C and 40°C (50°F and 104°F).
The display shows "Er5"	The thermometer is not working properly.	Take the batteries out and wait one minute. Insert the batteries and take a new measurement. If the display shows "Er5" again contact your retailer.

N.B. Please Note: When the device does not recover proper operation despite having performed these checks, contact the retailer.

DISPOSAL PROCEDURE (Dir. 2012/19/Eu-WEEE)

The symbol on the bottom of the device indicates the separated collection of electric and electronic equipment. At the end of life of the device, do not dispose it as mixed solid municipal waste, but dispose it referring to a specific collection centre located in your area or returning it to the distributor, when buying a new device of the same type to be used with the same functions.

If the appliance to be disposed of is less than 25 cm, it can be returned to a retail location that is over 400 m² without having to purchase a new, similar device. This procedure of separated collection of electric and electronic devices is carried out forecasting a European environmental policy aiming at safeguarding, protecting and improving environment quality, as well as avoiding potential effects on human health due to the presence of hazardous substances in such equipment or to an improper use of the same or of parts of the same.

Caution! The wrong disposal of electric and electronic equipment may involve sanctions. To correctly dispose of batteries (Dir. 2013/56/Eu) do not throw them into household waste but rather dispose of them as special waste at specialised recycling centres. For more information about disposal of used batteries, contact the store where you purchased the equipment containing the batteries, the municipality or the local waste disposal service.

WARRANTY

This equipment is covered by a 2-year warranty starting from the purchase date, which has to be proven by means **of the stamp or signature of the dealer and by the fiscal payment receipt that shall be kept enclosed to this document.** This warranty period complies with the laws in force and applies only in case the consumer is an individual. Laica products are

designed for domestic use only; therefore, their use for commercial purposes is forbidden. The warranty covers only production defects and does not apply, if the damage is caused by accidental events, misuse, negligence or improper use of the product.

Use the supplied accessories only; the use of other accessories invalidate the warranty cover. Do not open the equipment for any reason whatsoever; the opening or the tampering of the equipment, definitely invalidate the warranty. Warranty does not apply to parts subject to tear and wear along the use, as well as to batteries, when they are included. Warranty decays after 2 years from the purchase date. In this case, technical service interventions will be carried out against payment. Information on technical service interventions, under warranty terms or against payment, may be requested by writing at info@laica.com. Repair or replacement of products covered by the warranty terms are free of charge. In case of faults, refer to the dealer; DO NOT send anything directly to LAICA. All interventions covered by warranty (including product replacement or the replacement of any of its parts) shall not be meant to extend the duration of the original warranty period of the replaced product. The manufacturer cannot be made liable for any damages that may be caused, directly or indirectly, to people, things and pets resulting from the inobservance of the provisions contained in the suitable user manual, especially with reference to the instructions related to installation, use and maintenance of the equipment. Being constantly involved in improving its own products, Laica reserves the right to modify its product, totally or partially, with reference to any production need, without that this involves any responsibility whatsoever by Laica or its dealers.

STANDARDS

The device complies with the following standards: ISO80601-2-56 Medical electrical equipment - Part 2-56: Particular requirements for basic safety and essential performance of clinical thermometers for body temperature

measurement, EN 60601-1-11 Medical electrical equipment - Part 1-11: General requirements for basic safety and essential performance - Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment and complies with the requirements of EN 60601-1-2(EMC), IEC/EN60601-1(Safety) standards.

Electromagnetic compatibility
Table 1

Manufacturer's declaration and guidelines electromagnetic emissions		
The Device is intended for use in the electromagnetic environments specified below. The customer or user of the Device must make sure that the device is used in an appropriate environment.		
Emissions test	Compliance	Electromagnetic environment - guide
Radio frequency emissions CRISP 11	Group 1	The device only uses radio frequency energy for internal operation. Radio frequency emissions are therefore very low and do not cause any interference in electronic equipment located nearby.
Radio frequency emissions CRISP 11	Class B	The device is suitable for use in all domestic buildings and in those directly connected to the public low-voltage electricity supply network, intended for residential homes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage variations/flicker emissions IEC 61000-3-3	Not applicable	

Attention:

- This equipment should not be used adjacent to or stacked with other equipment as this may lead to malfunction. If such use is necessary, this equipment and the other equipment must be examined to check that it works normally.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) must not be used at a distance of less than 30 cm (12 inches) from any part of the ME ME EQUIPMENT or ME SYSTEM, including the cables specified by the producer. Otherwise the equipment may experience a decline in performance.

Table 2

Manufacturer's declaration and guidelines electromagnetic emissions			
The Device is intended for use in the electromagnetic environments specified below. The customer or user of the Device must make sure that the device is used in an appropriate environment.			
Immunity test	ICE 60601 test level	Level of compliance	Electromagnetic environment - guide
Electrostatic discharges (EDS) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors should be wood, concrete or ceramic. If the floors are covered with synthetic material, the relative humidity must be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for supply lines ±1 kV for input/output lines	Not applicable	
Overvoltage IEC 61000-4-5	± 1 kV line to line ± 2 kV line to ground	Not applicable	
Voltage dips, short interruptions and voltage variations on the power input lines IEC 61000-4-11	<5% UT (>95% dip in UT) for 0.5 cycles 40% UT (>60% dip in UT) for 5 cycles 70% UT (>60% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 seconds	Not applicable	
Network frequency magnetic field (50 Hz/60 Hz) IEC 61000-4-8	30 A/m; 50 Hz or 60 Hz	30 A/m; 50 Hz or 60 Hz	Network frequency magnetic fields should have characteristic levels of a typical location in a commercial or hospital environment.

Table 3

Manufacturer's declaration and guidelines – electromagnetic Immunity			
The Device is intended for use in the electromagnetic environments specified below. The customer or user of the Device must make sure that the device is used in an appropriate environment.			
Immunity test	ICE 60601 test level	Level of compliance	Electromagnetic environment - guide
Conducted radiofrequency IEC 61000-4-6	3 V from 0.15 MHz to 80 MHz 6 V from 0.15 MHz to 80 MHz outside ISM ranges	Not applicable	Portable and mobile radiofrequency communication devices must not be used at a distance "d" of any part of the Device, including cables, lower than the separation distance calculated from the equation relevant to the transmitter frequency. Recommended separation distance $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$ 80 MHz to 800 MHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$ 800 MHz to 2,7 GHz Where P is the maximum transmitter output index in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed radiofrequency transmitters, as determined by an electromagnetic site survey ^a , must be lower than the compliance level in each frequency range b. Interference may occur near equipment marked with the following symbol: 
Conducted radiofrequency IEC 61000-4-3	10 V/m From 80 MHz to 2.7 GHz	10 V/m	
Wireless radiofrequency IEC 61000-4-3	380 MHz, 27v/m 450 MHz, 28V/m 710 MHz, 745 MHz, 780 MHz 9 V/m 810 MHz, 870 MHz, 930 MHz 28V/m 1720 MHz, 1845 MHz, 1970 MHz 28V/m 5240 MHz, 5500 Mhz, 5785 MHz 9V/m	380 MHz, 27v/m 450 MHz, 28V/m 710 MHz, 745 MHz, 780 MHz 9 V/m 810 MHz, 870 MHz, 930 MHz 28V/m 1720 MHz, 1845 MHz, 1970 MHz 28V/m 5240 MHz, 5500 Mhz, 5785 MHz 9V/m	

Table 4

Separation distances recommended between portable and mobile RF communication devices and the Device.		
The device is intended for use in electromagnetic environments in which radiated radiofrequency disturbances are controlled. The customer or the user of the Device can contribute to the prevention of electromagnetic interference by maintaining a minimum distance between the portable and mobile radiofrequency communication devices (transmitters) and the Device based on the indications provided below, according to the maximum output power of communication devices.		
Maximum rated output power of the transmitter (W)	Separation distance depending on the frequency of the transmitter (m)	
	From 80 kHz to 800 MHz $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$	From 800 MHz to 2,7 GHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0,01	0,12	0,23
0,1	0,38	0,73
1	1,2	2,3
10	3,8	7,3
100	12	23
For transmitters whose maximum rated output power is not among those listed above, the recommended separation distance d in meters (m) can be determined by the equation relevant to the transmitter frequency, where P is the maximum rated output power of the transmitter, expressed in watts (W) according to the data provided by the transmitter manufacturer. NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is influenced by absorption and reflection of structures, objects and people.		

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NORMAS

El dispositivo respeta las normas: ISO80601-2-56 Medical electrical equipment - Part 2-56:Particular requirements for basic safety and essential performance of clinical thermometers for body temperature measurement, EN 60601-1-11 Medical electrical equipment - Part 1-11: General requirements for basic safety and essential performance - Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment and complies with the requirements of EN 60601-1-2(EMC), IEC/EN60601-1(Safety) standards.

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NORMAS

O dispositivo está em conformidade com as normas: ISO80601-2-56 Medical electrical equipment - Part 2-56:Particular requirements for basic safety and essential performance of clinical thermometers for body temperature measurement, EN 60601-1-11 Medical electrical equipment - Part 1-11: General requirements for basic safety and essential performance - Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment and complies with the requirements of EN 60601-1-2(EMC), IEC/EN60601-1(Safety) standards.

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NORMEN

Die Vorrichtung entspricht den Normen: ISO80601-2-56 Medical electrical equipment - Part 2-56:Particular requirements for basic safety and essential performance of clinical thermometers for body temperature measurement, EN 60601-1-11 Medical electrical equipment - Part 1-11: General requirements for basic safety and essential performance - Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment and complies with the requirements of EN 60601-1-2(EMC), IEC/EN60601-1(Safety) standards.

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NORMES

Le dispositif est conforme aux normes: ISO80601-2-56 Medical electrical equipment - Part 2-56:Particular requirements for basic safety and essential performance of clinical thermometers for body temperature measurement, EN 60601-1-11 Medical electrical equipment - Part 1-11: General requirements for basic safety and essential performance - Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment and complies with the requirements of EN 60601-1-2(EMC), IEC/EN60601-1(Safety) standards.

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Η συσκευή είναι σύμφωνη με τους κανονισμούς: ISO80601-2-56 Medical electrical equipment - Part 2-56:Particular requirements for basic safety and essential performance of clinical thermometers for body temperature measurement, EN 60601-1-11 Medical electrical equipment - Part 1-11: General requirements for basic safety and essential performance - Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment and complies with the requirements of EN 60601-1-2(EMC), IEC/EN60601-1(Safety) standards.

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STANDARDE

Dispozitivul corespunde normelor: ISO80601-2-56 Medical electrical equipment - Part 2-56:Particular requirements for basic safety and essential performance of clinical thermometers for body temperature measurement, EN 60601-1-11 Medical electrical equipment - Part 1-11: General requirements for basic safety and essential performance - Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment and complies with the requirements of EN 60601-1-2(EMC), IEC/EN60601-1(Safety) standards.

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NORMY

Tento produkt splňuje následující normy: ISO80601-2-56 Medical electrical equipment - Part 2-56:Particular requirements for basic safety and essential performance of clinical thermometers for body temperature measurement, EN 60601-1-11 Medical electrical equipment - Part 1-11: General requirements for basic safety and essential performance – Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment and complies with the requirements of EN 60601-1-2(EMC), IEC/EN60601-1(Safety) standards.

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STANDARDE

Izdelek izpolnjuje naslednje standarde: ISO80601-2-56 Medical electrical equipment - Part 2-56:Particular requirements for basic safety and essential performance of clinical thermometers for body temperature measurement, EN 60601-1-11 Medical electrical equipment - Part 1-11: General requirements for basic safety and essential performance – Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment and complies with the requirements of EN 60601-1-2(EMC), IEC/ EN60601-1(Safety) standards.

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SZABVÁNYOK

A keszulek az alábbi szabványoknak felel meg: ISO80601-2-56 Medical electrical equipment - Part 2-56:Particular requirements for basic safety and essential performance of clinical thermometers for body temperature measurement, EN 60601-1-11 Medical electrical equipment - Part 1-11: General requirements for basic safety and essential performance – Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment and complies with the requirements of EN 60601-1-2(EMC), IEC/EN60601-1(Safety) standards.

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STANDARDI

Izdelek ustreza naslednjim standardom: ISO80601-2-56 Medical electrical equipment - Part 2-56:Particular requirements for basic safety and essential performance of clinical thermometers for body temperature measurement, EN 60601-1-11 Medical electrical equipment - Part 1-11: General requirements for basic safety and essential performance – Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment and complies with the requirements of EN 60601-1-2(EMC), IEC/EN60601-1(Safety) standards.

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Za dodatne informacije: www.laica.it

STANDARDI

Proizvod ispunjava sljedeće standarde: ISO80601-2-56 Medical electrical equipment - Part 2-56:Particular requirements for basic safety and essential performance of clinical thermometers for body temperature measurement, EN 60601-1-11 Medical electrical equipment - Part 1-11: General requirements for basic safety and essential performance – Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment and complies with the requirements of EN 60601-1-2(EMC), IEC/EN60601-1(Safety) standards.

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