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GB	English	FACIAL SAUNA Instructions and warranty
-----------	----------------	--

YOU ARE STRONGLY RECOMMENDED TO READ THE INSTRUCTIONS AND WARNINGS HEREIN BEFORE USING THE EQUIPMENT AND TO KEEP THE BOOKLET IN A SAFE PLACE.

The skin reflects the overall health status of a person. It is possible to try to improve it following a balanced diet, performing a regular physical activity, and sleeping well in the night.

Stress and nervous tension, digestion problems, weakening of the immune system, and bad blood circulation can emphasize skin problems like dermatitis, acne, herpes, and psoriasis.

It is important to treat one's own skin very delicately, avoiding the use of irritating substances. In the morning, wash face and neck with a delicate detergent. Rinse carefully. Use a tonic without alcohol to close the pores and then, dry it slightly. Then apply a moistened cotton pad. Let it dry and then apply some moisturizing cream. Before performing these operations, place the mirror in the most comfortable position for you to carry them out more easily. Apply masks specific for your type of skin, peeling and exfoliating products at least once a week.

The face sauna by LAICA, exploiting the beneficial action of the steam, allows a deep and accurate cleaning of the face removing any impurities and blackheads. The refreshed and moisturized skin will seer younger and brighter. Ideal for use in the home, it removes makeup deep down, opening the skin's pores and laying the foundations for a good tan as it helps prevent the skin going dry. Comes complete with measuring beaker.

The face sauna by LAICA, exploiting the beneficial action of the steam, allows a deep and accurate cleaning of the face removing any impurities and blackheads. The refreshed and moisturized skin will seer younger and brighter. Ideal for use in the home, it removes makeup deep down, opening the skin's pores and laying the foundations for a good tan as it helps prevent the skin going dry.

Comes complete with measuring beaker.

GENERAL WARNINGS

Using any electrical appliance entails compliance with a number of fundamental rules:

- Unplug the appliance from the mains straight after use and before commencing cleaning or inserting/removing accessories.
- Do not unplug by tugging on either power cord or appliance.
- If an electrical appliance falls into the water, do not attempt to retrieve it; instead, unplug immediately.
- Do not use the appliance while bathing or showering.
- Never immerse the appliance in water or other liquids.
- Children, invalids or the disabled must not use the appliance when unattended. Do not allow children to play with the device.
- This appliance must be used solely for the purpose for which it was designed and in accordance with the operating instructions. Any other use may cause damage to the appliance and hence hazardous. The manufacturer cannot be held responsible for any damage resulting from improper or incorrect use.
- Before using the appliance, make sure it is intact and does not feature any visible signs of damage. If in any doubt, do not use the appliance and contact your dealer.
- In the event of a fault and/or irregular operation, turn the appliance off at the main supply and contact your dealer for repairs.
- Keep the appliance away from sources of heat.
- Make sure your hands are dry when unplugging or plugging in the appliance and when turning the function control knob.
- This appliance must only be used with the voltage given on the bottom of the product.
- This appliance is not suitable for professional or medical use.

After using the appliance for 15 minutes non-stop, it must be allowed to cool for at least 15 minutes.

- Wait for few minutes before the facial sauna starts to give off the steam.
- Do not keep your face in the same position too near the steam diffuser for too long as this might lead to scalding.
- Do not use the device close to water containers, showers, and wash-basins.
- Should the power supply cable be damaged, it is necessary that the manufacturer repairs it. In any case, for repairings always contact your dealer.

IMPORTANT

- Turn the apparatus off immediately and consult a doctor should the treatment be unpleasant or cause pain.

ATTENTION:

Always consult your doctor in case of:

- **fever**
- **purulent diseases**
- **inflammation**
- **pain**
- **if using medicines or anaesthetics**
- **skin lesions or diseases affecting the face.**

HOW TO USE THE FACIAL SAUNA

- A Wash your face well before using the apparatus.
- B Place the facial sauna on a flat, hard surface.
- C Pour warm water in the small basin of the main unit using the proper measuring device and making sure not to exceed the indicated maximum level (MAX). Warning: If the water accidentally overflows out of the tray in the main unit, empty everything, dry it carefully and fill it up again.
- D After selecting the required accessory (facial cone or inhaler), connect it to the main unit.
- E Plug the appliance into the wall socket (230V).
- F Switch on the appliance by setting the switch to "1". The orange light indicates the appliance is on; wait 3/4 minutes for the water to heat.
- G Close your eyes and relax.
- H When all the water in the tank has evaporated, switch off the appliance and allow it to cool for at least 15 minutes.
- I Switch off the appliance by setting the switch to "0".
- J Remove the steam diffuser or inhaler from the main unit.
- K Unplug the appliance from the mains and empty any water left in the tank.

INHALER

In addition to the steam diffuser used for inhaling through the mouth, there is a special mask for inhaling through the nose. When it is cold, mucous membranes become particularly dry. Using the appliance with water, the membranes are moistened so that they work better, i.e., purify the air we breathe. If symptoms persist, consult a doctor. Adding additives to the water might have various effects. Act on advice given by your doctor or chemist only.

ACCESSORIES/COMPONENTS

- 1) Steam diffuser
- 2) Water tank
- 3) Main unit
- 4) Switch
- 5) Lamp indicator light
- 6) Electric cable
- 7) Measuring beaker
- 8) Inhaler

MAINTENANCE

Wait until the apparatus cools down before carrying out any maintenance.

- Clean the small water basin and face mask with a damp cloth after each use to avoid the deposit of limestone or other residuals.
- Clean the main unit with a dry cloth.
- Do not use chemicals or abrasive products for cleaning.
- Do not immerse the main unit in water or other liquids.

GB	English	FACIAL SAUNA Instructions and warranty
-----------	----------------	--

This product is meant for domestic use.

The conformity, certified by the CE-marking on the device, refers to the EEC Directive 2004/108 on electromagnetic compatibility and to the EEC Directive 2006/95 on electric material destined to be used within given voltage limits. This product is an electronic device that was tested to assure, according to the current technical knowledge, that it does not interfere with other nearby located equipment (electromagnetic compatibility) and that it is safe if used according to the advices given in the instructions for use. In case of operation anomalies, do not keep using the device and, if necessary, refer to the manufacturer directly.

DISPOSAL PROCEDURE (Dir. 2012/19/Eu-WEEE)

The symbol on the bottom of the device indicates the separated collection of electric and electronic equipment.

At the end of life of the device, do not dispose it as mixed solid municipal waste, but dispose it referring to a specific collection centre located in your area or returning it to the distributor with your purchase. Wash it using an exfoliating, yet delicate, product and then, dry it slightly. Then apply a moistened cotton pad. Let it dry and then apply some moisturizing cream. Before performing these operations, place the mirror in the most comfortable position for you to carry them out more easily. Apply masks specific for your type of skin, peeling and exfoliating products at least once a week. The face sauna by LAICA, exploiting the beneficial action of the steam, allows a deep and accurate cleaning of the face removing any impurities and blackheads. The refreshed and moisturized skin will seer younger and brighter. Ideal for use in the home, it removes makeup deep down, opening the skin's pores and laying the foundations for a good tan as it helps prevent the skin going dry.

WARRANTY

This equipment is covered by a 2-year warranty starting from the purchase date, which has to be proven by means of the **stamp or signature of the dealer and by the fiscal payment receipt that shall be kept enclosed to this document.** This warranty period complies with the laws in force and applies only in case the consumer is an individual. Laica products are designed for domestic use only; therefore, their use for commercial purposes is forbidden.

The warranty covers only production defects and does not apply, if the damage is caused by accidental events, misuse, negligence or improper use. The manufacturer cannot be held liable for any damages that may be caused, directly or indirectly, to people, things and pets resulting from the inobservance of the provisions contained in the suitable user manual, especially with reference to the instructions related to installation, use and maintenance of the equipment.

Being constantly involved in improving its own products, Laica reserves the right to modify its product, totally or partially, with reference to any production or model without that this involves any responsibility whatsoever by Laica or its dealers.

Laica products are designed for domestic use only; therefore, their use for commercial purposes is forbidden. The warranty covers only production defects and does not apply, if the damage is caused by accidental events, misuse, negligence or improper use. The manufacturer cannot be held liable for any damages that may be caused, directly or indirectly, to people, things and pets resulting from the inobservance of the provisions contained in the suitable user manual, especially with reference to the instructions related to installation, use and maintenance of the equipment. Being constantly involved in improving its own products, Laica reserves the right to modify its product, totally or partially, with reference to any production or model without that this involves any responsibility whatsoever by La