

	<b>BM7002</b> 	<b>BM7003</b> 	<b>BM2605</b> 	<b>BM2604</b> 
Tipologia Type	Misuratore da braccio Upper arm BPM	Misuratore da polso Wrist BPM	Misuratore da braccio Upper arm BPM	Misuratore da braccio Upper arm BPM
Circonferenza braccio/polso Arm/wrist circumference	standard 220 - 300 mm	standard 140 - 195 mm	extra-large 220 - 420 mm	standard 220 - 300 mm
Smart	●	●		
Capacità memorie Memory capacity	(60 memories x 2 users)	(60 memories x 2 users)	90 memories	120 (30 memo- ries x 4 users)
Media delle misurazioni Results average	●	●	●	●
Data e ora Date and time	●	●	●	●
Rilevazione battiti irregolari Irregular heartbeat detection	●	●	●	●
Confronto dati con valori OMS Comparison of data to WHO values	●	●	●	●
Display	LCD - 6.3 X 8 cm	LCD - 4.2 x 4.7 cm	LCD - 8 x 7 cm	LCD - 8.7 x 7.4 cm
Custodia Storage case	●	●	●	●
Alimentazione Power supply	4 x 1.5V alkaline batteries AA (included)	2 x 1.5V alkaline batteries AAA (included)	4 x 1.5V alkaline batteries AA (included)	4 x 1.5V alkaline batteries AAA (included)
Dimensioni Size	15 x 9 x 5.1 cm	6.2 x 8.1 x 6.4 cm	16.5 x 23.5 x 10	15 x 9.5 x 4.5 cm
Peso Weight	472 g	190 g	750 g	420 g

	<b>BM2302</b> 	<b>BM2301</b> 	<b>BM2006</b> 	<b>BM1003</b> 	<b>BM1007</b> 	<b>BM1006</b> 
Tipologia Type	Misuratore da braccio Upper arm BPM	Misuratore da braccio Upper arm BPM	Misuratore da braccio Upper arm BPM	Misuratore da polso Wrist BPM	Misuratore da polso Wrist BPM	Misuratore da polso Wrist BPM
Circonferenza braccio/polso Arm/wrist circumference	standard 220 - 300 mm	standard 220 - 360 mm	standard 220 - 360 mm	standard 140 - 195 mm	standard 140 - 195 mm	standard 135 - 215 mm
Smart						
Capacità memorie Memory capacity	120 (30 memories x 4 users)	120 (30 memories x 4 users)	120 memories	60 memories	120 (30 memories x 4 users)	(60 memories x 2 users)
Media delle misurazioni Results average	●	●			●	●
Data e ora Date and time	●	●			●	●
Rilevazione battiti irregolari Irregular heartbeat detection	●	●			●	●
Confronto dati con valori OMS Comparison of data to WHO values	●	●			●	●
Display	LCD - 5.1 x 6.1 cm	LCD - 5.5 x 8.2 cm	LCD - 4.5 x 3 cm	LCD backlit 4.6 x 6.8 cm	LCD - 4 x 4.7 cm	LCD - 3.7 x 4.8 cm
Custodia Storage case	●	●	●	●	●	●
Alimentazione Power supply	4 x 1.5V alkaline batteries AAA (included)	4 x 1.5V alkaline batteries AA (included)	4 x 1.5V alkaline batteries AA (included)	2 x 1.5V alkaline batteries AAA (included)	2 x 1.5V alkaline batteries AAA (included)	2 x 1.5V alkaline batteries AAA (included)
Dimensioni Size	10.6 x 8.2 x 5.2 cm	16.3 x 11 x 8.2 cm	13.2 x 10 x 6.5 cm	8.8 x 7.9 x 6.8 cm	6.8 x 9 x 7 cm	7.7 x 6.2 x 6.5 cm
Peso Weight	380 g	530 g	450 g	260 g	190 g	180 g